

INTRODUCING FORTRESS DIAGNOSTICS GEL CARDS

“An innovative method of testing”

NEW

Key Benefits:



Reduced working time

The gel card method vastly reduces preparation time due to multi testing and no requirement for washings steps.



Cost Effective

Fortress Diagnostics supply gel cards at a competitive price promoting this method as an innovative solution in preference to conventional methods.



Accurate Results



Fortress Diagnostics gel cards are proven to provide accurate results in minimal time.



Complete Solution

One test provides comprehensive results for blood grouping compared to multiple tests in conventional tube methods.

New lines coming soon!



Fortress
Diagnostics

www.fortressdiagnostics.com



DISCOVER MORE

Fortress Forward Group 1



Identifies ABO and Rh D groups by using gel method. One card involves 2 tests. Includes control well.

PACKING	CAT NO.
12 Fortress Forward Group 1 Gel Cards	BGGCFG01-12
48 Fortress Forward Group 1 Gel Cards	BGGCFG01-48

Fortress Forward Group 2



Identifies ABO and Rh D groups by using gel method. One card involves 2 tests. Does not include control well, instead includes AB.

PACKING	CAT NO.
12 Fortress Forward Group 2 Gel Cards	BGGCFG02-12
48 Fortress Forward Group 2 Gel Cards	BGGCFG02-48

Fortress Forward Group 3



Provides the detection of subgroups of A and variants of D on the same gel card.
Identifies DVI- weak reactive and partial variants of D but not the VI variant of D directly.
Identifies DVI+ reactive, variant of DVI includes D variants and weak D's. Identifies the erythrocytes of antibody Rh DVI directly.

PACKING	CAT NO.
12 Fortress Forward Group 3 Gel Cards	BGGCFG03-12
48 Fortress Forward Group 3 Gel Cards	BGGCFG03-48

Fortress Forward and Reverse Grouping



Provides the detection of ABO, Rh D, Rh D VI ve reverse grouping on one gel card.
Identifies DVI- weak reactive and partial variants of D but not the VI variant of D directly.
Identifies DVI+ reactive, variant of DVI includes D variants and weak D's. Identifies the erythrocytes of antibody Rh DVI directly.
At the last two wells A1 and B cells are used for reverse apply.

PACKING	CAT NO.
12 Fortress Forward and Reverse Grouping Gel Cards	BGGCFRG1-12
48 Fortress Forward and Reverse Grouping Gel Cards	BGGCFRG1-48


Fortress Comp Test



The gel card for forward grouping and Cross-match
Provides detection of ABO, Rh mismatch, the existence of antibodies that are against donators erythrocytes in patients serum or the patients erythrocytes in donators serum.

PACKING	CAT NO.
12 Fortress Comp Test Gel Cards	BGGCFCT1-12
48 Fortress Comp Test Gel Cards	BGGCFCT1-48

Fortress Neonatal Grouping



Provides applying of ABO, RH D, Rh D VI and direct Coombs test on newborns on one gel card.
AHG well involves Polyclonal anti-IgG and monoclonal C3 antibodies, IgG well involves monospecific Anti-IgG antibodies.

PACKING	CAT NO.
12 Fortress Neonatal Grouping Gel Cards	BGGCNEOG-12
48 Fortress Neonatal Grouping Gel Cards	BGGCNEOG-48


Fortress Rh Phenotype +K



Provides identifying of Rh and K antibodies which have clinic importance by gel method.

PACKING	CAT NO.
12 Fortress Rh Phenotype +K Gel Cards	BGGCRHPK-12
48 Fortress Rh Phenotype +K Gel Cards	BGGCRHPK-48

Fortress A subgroup



This test provides identifying of A1 antibody inside erythrocytes by using gel method.

PACKING	CAT NO.
12 Fortress A Subgroup Gel Cards	BGGCASUB-12
48 Fortress A Subgroup Gel Cards	BGGCASUB-48

Fortress Neutral Gel



Provides applying of reverse grouping, antibody screening and determining, compability tests.

PACKING	CAT NO.
12 Fortress Neutral Gel Cards	BGGCNGEL-12
48 Fortress Neutral Gel Cards	BGGCNGEL-48

Fortress AHG Gel



Provides applying of Coombs test (DAT) directly, Coombs test (IAT), compatibility, antibody screening and determining tests indirectly. All wells include Polyspecific anti-human globulin.

PACKING	CAT NO.
12 Fortress AHG Gel Cards	BGGCAHGG-12
48 Fortress AHG Gel Cards	BGGCAHGG-48

Fortress DAT Screen



Provides determining the erythrocytes that get sensitive as in-vivo with IgG or complement C3d.
Involves IgG well monospecific Anti-IgG, C3d well Anti human C3, AHG well Polyspecific anti-human globuin.

PACKING	CAT NO.
12 Fortress DAT Screen Gel Cards	BGGCDATS-12
48 Fortress DAT Screen Gel Cards	BGGCDATS-48

Fortress Neutral/AHG Gel



Provides applying of Indirect Coombs test (IAT), compatibility, antibody screening and determining tests indirectly.

PACKING	CAT NO.
12 Fortress Neutral/AHG Gel Cards	BGGCNAHG-12
48 Fortress Neutral/AHG Gel Cards	BGGCNAHG-48

OTHER ITEMS	SIZE	CAT NO.
Fortress Bromelain Solution	10ml	BGGCBS10
Fortress Bromelain Solution	50ml	BGGCBS50
Fortress LISS Diluent	100ml	BGGCL100
Fortress LISS Diluent	500ml	BGGCL500
Fortress Papain Solution	10ml	BGGCPAPS
Fortress Cell Reverse A1, B	-	BGGCRBC2
Fortress Cell Reverse A1, A2, B, O	-	BGGCRBC4
Fortress Cell Screen 1-2	-	BGGCSC12
Fortress Cell Screen 1-2 Pap	-	BGGCSC2P
Fortress Cell Screen 1-2-3	-	BGGCSC23