

# GEL CARD CENTRIFUGE

CAT. NO.: CENF0024

## Key Benefits:



### Compact Solution

Benchtop Centrifuge with a max speed of 3,000 RPM. Only 58x17x38 CMS.



### Easy to Use

This Centrifuge is very easy to use and has a built in LCD display with 2 lines of text with 16 characters.



### Time Control

Set time and walk away with two progressive timers < 99 min or continuous setting.




### Cost Effective

Low startup cost for small to medium labs who require a centrifuge for their Gel Cards. Fortress Diagnostics also supply gel cards at a competitive price.

**Reader and Incubator also available!**



[www.fortressdiagnostics.com](http://www.fortressdiagnostics.com)



**Fortress**  
Diagnostics



DISCOVER MORE

## TECHNICAL SPECIFICATIONS

DESCRIPTION	SPECIFICATION
Power Supply:	220 V / 50 – 60 Hz
Power requirement:	100 W
Max. rotational speed:	3,000 rpm
Max. Load:	24 Gel Card
Display:	2 line x 16 character LCD Blue
Time control:	Two Progressive, timed < 99 min or continuous
Time counting modes:	Selectable, at set speed or from starting
Ambient temperature:	2 °C – 35 °C
Max. relative humidity:	75 %
Safety lid lock:	Yes
Lid drop protection:	Yes
Opening the Lid:	Manual
Dimensions:	Height: 580 mm / Depth: 170 mm / Width: 380 mm
Weight (without rotor):	16.0 kg
Acceleration time (230 V):	< 13 sec
Deceleration time:	< 11 sec.
Max Noise (dBA):	<65 dBA

